

**Preparation and Titration of HIV-1 Env-Pseudoviruses  
(Montefiori Lab)  
July 2007**

## **I. INTRODUCTION**

The performance of neutralizing antibody assays under properly standardized, optimized and validated conditions requires accurate titration of virus infectivity. Moreover, use of molecularly cloned pseudoviruses has advantages over uncloned virus for greater reagent stability and assay reproducibility and precision. This protocol describes the production of molecularly cloned Env-pseudotyped HIV1 in 293T/17 cells by co-transfection with an Env-expressing plasmid plus a backbone plasmid lacking Env. Co-transfection generates pseudovirus particles that are able to infect cells but, due to the absence of a complete genome, are unable to produce infectious progeny virions. This single round of infection is readily detectable in genetically engineered cells lines that contain a Tat-responsive reporter gene, such as luciferase. Neutralization assays based on a single-round of infection should use viruses that are titrated in a similar single-round infection format.

## **II. DEFINITIONS**

GM, Growth Medium

DMEM, Dulbecco's Modified Eagle Medium

HEPES, N-2-Hydroxyethylpiperazine-N'-2-Ethanesulfonic Acid

Luc, Luciferase

RLU, Relative Luminescence Units

FBS, Fetal Bovine Serum

ID, Identification

## **III. REAGENTS AND MATERIALS**

Recommended vendors are listed. Unless otherwise specified, products of equal or better quality may be used when necessary.

### **293T/17 cells**

ATCC

### **TZM-bl cells**

NIH AIDS Research and Reference Reagent Program

### **Growth Medium\***

**DMEM, with L-glutamine, sodium pyruvate, glucose and pyridoxine, sterile. Store refrigerated at 4°C.**

Gibco BRL Life Technologies, cat no. 11995-065

**Fetal bovine serum, heat-inactivated 56°C for 30 minutes, 500 ml bottle, sterile. Store at -20°C. Once thawed, store at 4°C for up to 1 month.**

Hyclone

### **Gentamicin solution, 10 mg/ml, sterile**

Sigma, store at 4°C

### **HEPES Buffer, 1M, store at 4°C**

Sigma

\*Complete GM consists of DMEM containing 10% heat-inactivated FBS, 50 µg gentamicin/ml and 25mM HEPES. To make 500 ml of GM, combine 435 ml DMEM, 50 ml FBS, 2.5 ml of gentamicin, and 12.5 ml HEPES in a sterile bottle, mix, store at 4°C for up to 2 months. Warm medium to 25°C – 37°C prior to use.

**DEAE dextran, hydrochloride, average Mol. Wt. 500,000**

Sigma

Prepare a 7.5 mg/ml solution by dissolving 3.75 gm in 500 ml of sterile water. Store at -80°C in 10 ml aliquots in 15 ml sterile polypropylene tubes.

**Trypsin-EDTA (0.25% trypsin, 1 mM EDTA), sterile**

Gibco cat. no. 25200-056

**FuGENE 6 Transfection Reagent**

Roche Applied Science cat no. 11988387001

**Britelite Luminescence Reporter Gene Assay System**

Perkin Elmer Life and Analytical Sciences cat no. 6016979

Reconstitute one vial of lyophilized Britelite Substrate Solution with 250 ml of Britelite Substrate Buffer Solution. After the substrate has dissolved completely (about 1 minute), mix gently and distribute 10.5 ml to 15 ml conical polypropylene tubes and store at -80°C immediately. Thaw in a room temperature water bath in the dark immediately before each use. Mix gently prior to use. Use within 60 minutes of thawing. Excess reagent may be stored at -80°C and used once more.

*Caution: The lyophilized Britelite substrate is classified as hazardous. Latex gloves, surgical gown and eye protection are required when working with these reagents.*

**Hemocytometer**

Hausser Scientific, Horsham, PA

**12-channel pipettman, 5-50 µl**

ThermoLabsystem

**12-channel pipettman, 30-300 µl**

ThermoLabsystem

**Single channel pipettman, 5-50 µl**

ThermoLabsystem

**Single channel pipettman, 30-200 µl**

ThermoLabsystem

**Microliter pipet tips, sterile**

ICN

**PipetteAid XP**

Drummond Scientific Co.

**Disposable pipettes, sterile, individually wrapped**

Falcon/VWR

1 ml pipets

5 ml pipets

10 ml pipets

25 ml pipets

50 ml pipets

**Flat-bottom culture plates, 96-well, low evaporation, sterile**

Costar/VWR

**Flat-bottom black solid plates, 96-well, Costar brand**

Fisher

**15 ml conical polypropylene tubes, screw-cap, sterile**

Corning

**Culture flasks with vented caps, sterile**

Costar/VWR

T-25 flask

T-75 flask

**Reagent reservoirs, 50 ml capacity**

Costar

**Instrumentation:**

**Biological Laminar Flow Cabinet (annual certification required)**

NuAIRE, Plymouth, MN; Model NU-425-600

**Incubator, water-jacketed (37°C, 5% CO<sub>2</sub> standard requirements)**

Forma Scientific, Steri-Cult 200, Model 3033

**Desk-top centrifuge (low speed capable of up to 500 x g)**

Jouan, Model C412

50 ml tube holder

15 ml tube holder

Microtitration plate holder

**Luminometer equipped to read 96-well plates**

PerkinElmer Life Sciences, Model Victor2

**Controlled temperature water bath**

Precision Scientific, Model 182

**Microcentrifuge (maximum rotational speed = 14,000 rpm)**

Eppendorf, Model 5415C

18 place standard rotor F-45-18-11 for 1.5 ml microcentrifuge tubes

**Specimens:**

Molecularly cloned pseudoviruses are generated by using a two plasmid system: 1) an Env expression plasmid (e.g., pcDNA 3.1D/V5-His-TOPO-Env) and a backbone vector (e.g., pSG3ΔEnv) that expresses the entire HIV-1 genome except Env. Unclassified viruses for titration may be cell-free stocks that are produced in either PBMC or human T cell lines.

#### **IV. PROTOCOL**

*NOTE 1: All incubations are performed in a humidified 37°C, 5% CO<sub>2</sub> incubator unless otherwise specified.*

##### **Thawing Cells**

*NOTE 2: Be sure to wear a full-face shield during handling of frozen samples.*

*NOTE 3: TZM-bl and 293T/17 are adherent cell lines that are maintained in T-75 culture flasks. Cell monolayers are disrupted and removed by treatment with trypsin/EDTA at confluency.*

1. Transfer cryovials containing frozen cells from liquid nitrogen to a room temperature water bath in the biosafety hood. If liquid nitrogen has seeped into the cryovial, loosen the cap slightly to allow the nitrogen to escape during thawing. Hold the cryovial on the surface of the water bath with an occasional gentle “flick” during thawing. Do not leave the cryovial unattended during the thawing process. (It is important for cell viability that the cells are thawed and processed quickly - thawing only takes a few seconds). Dry off the outside of the cryovials and wipe with alcohol solution before opening to prevent contamination.
2. Transfer the contents of one vial of cells to a T-75 culture flask containing 30 ml of GM. Note: It is important to dilute the DMSO at least 30-fold at this point to avoid cell toxicity.
3. Incubate the cells at 37°C for 1 day.
4. Remove the medium and replace with 15 ml of fresh GM. Change the medium every 2-3 days until the cell monolayers are confluent.

##### **Splitting Cells**

1. Decant the culture medium and remove residual serum by rinsing monolayers with 5 ml of sterile PBS.
2. Slowly add 2.5 ml of an 0.25% Trypin-EDTA solution to cover the cell monolayer. Incubate at room temp for 30 seconds. Decant the trypsin solution and incubate at 37°C for 4 minutes in the case of TZM-bl cells. Incubate at room temperature for 1 minute in the case of 293T/17 cells. Do not agitate the cells by hitting or shaking the flask while waiting for the cells to detach.
3. Add 10 ml of GM and suspend the cells by gentle pipet action. Count cells.
4. Seed new T-75 culture flasks with approximately  $10^6$  cells in 15 ml of GM. Cultures are incubated at 37°C in a 5% CO<sub>2</sub>/95% air environment. Cells should be split approximately every 3 days.

##### **Transfection of 293T/17 Cells**

1. Seed  $3 \times 10^6$  293T/17 cells in a T75 flask containing 12 ml GM. Incubate overnight. Monolayers should be 50-80% confluent on the day of transfection.

2. Based on the volume of Env plasmid DNA and backbone plasmid DNA to be dispensed, add the appropriate volume of DMEM, such that the total volume of the mixture is 100  $\mu$ l, into one sterile tube. Dispense 4  $\mu$ g of Env plasmid DNA and 8  $\mu$ g of backbone plasmid DNA to the tube containing DMEM and mix well.
3. To a second sterile tube, add 652  $\mu$ l DMEM. Pipet 48  $\mu$ l of FuGENE 6 reagent directly into medium without contacting the sides of the plastic tube. Mix well.
4. Transfer the entire contents of the plasmid DNA mixture from the first tube to the second tube containing the FuGENE solution. Mix by pipetting or briefly vortexing.
5. Incubate for 30 minutes at room temperature to allow complex formation.
6. Add the entire contents of the transfection complexes to a T-75 flask of 293T/17 cells. Gently swirl the flask for uniform distribution of the complexes.
7. Incubate for 3 to 8 hours at 37°C in a 5% incubator to allow the plasmids to enter the cells.
8. Decant the medium containing DNA-FuGENE complexes and replace with 15 ml fresh GM. Incubate for 2 days.
9. Harvest virus-containing culture supernatants by collecting the medium from the flask using a pipet. Collect as much as possible without drawing cells into the pipet. Filter the virus-containing culture fluid through a 0.45-micron filter. Adjust the FBS concentration to 20%, mix. Distribute aliquots of appropriate volume to polypropylene screw-cap tubes that have been labeled to identify the isolate name and the date of harvest. The harvest date becomes the specific lot number. Store the aliquots at -80°C. Record the harvest and location of the vials. Include the identification regarding the Env plasmid DNA, backbone DNA, and cells used in the transfection when logging this information.
10. Add 12 ml of fresh GM to the cells in each flask. Incubate overnight and harvest the virus-containing culture supernatants once more, as previously indicated, and discard the cells. Record appropriate harvest information as before.

#### **Titration Virus in TZM-bl (TCID50 Assay)**

1. Place 100  $\mu$ l of GM per well in all wells of a 96-well flat-bottom culture plate. Transfer 25  $\mu$ l of virus to the first 4 wells of a dilution series (column 1, rows A-D for one virus and rows E-H for a second virus), mix, do serial 5-fold dilutions (i.e., transfer 25  $\mu$ l, mixing each time) for a total of 11 dilutions. Discard 25  $\mu$ l from the 11th dilution. Wells in column 12 will serve as cell controls (no virus added).
2. Add 100  $\mu$ l of TZM-bl cells (10,000 cells/100  $\mu$ l DMEM containing 30  $\mu$ g DEAE dextran/ml) to all wells. Rinse your pipet tips in a reservoir of RPMI between each step to minimize carry-over. The final concentration of DEAE dextran is 15  $\mu$ g/ml.

*NOTE 4: The concentrations of DEAE dextran shown above are approximations. The actual optimal concentration should be determined for each new batch of dextran by performing a titration assay: The dextran should be serially diluted in a 96-well plate containing GM, virus and TZM-bl cells added, and incubated for 48 hours. Following the incubation, luminescence should be measured and the appropriate concentration selected that yields the highest RLU but has no detrimental effect on the*

*cells as determined by microscopic examination.*

3. Incubate for 48 hours.

*Assays with replication-competent viruses must be incubated for 48 hours only to minimize virus replication.*

4. Remove 100  $\mu$ l of culture medium from each well, leaving approximately 100  $\mu$ l. Dispense 100  $\mu$ l of Britelite Reagent to each well. Incubate at room temperature for 2 minutes to allow complete cell lysis. Mix by pipet action (two strokes) and transfer 150  $\mu$ l to a corresponding 96-well black plate. Read the plate immediately in a luminometer.

5. Calculate the TCID<sub>50</sub> according to the method of Reed and Muench as described (8.1) using the "TCID<sub>50</sub>" macro. Wells with RLU <2.5 times background are considered negative for the calculation.

## **V. REFERENCES**

1. Johnson, V.A., and R.E. Byington. 1990. Infectivity assay (virus yield assay). In : Techniques in HIV Research (Aldovani, A., and Walker, B.D., eds.). Stockton Press, New York, N.Y., pp71-76.